

ACTIVITY REPORT

TITLE OF EXPERIMENT:

TITRATION OF HCL AGAINST NAOH

INTRODUCTION:


This is an experiment to titrate Hydrochloric Acid (HCl) against Sodium Hydroxide (NaOH) through the neutralisation reaction.

OBJECTIVE:

This experiment aims to discover what exact volume of an acid is needed to just neutralise a certain volume of an alkali.

MATERIALS AND APPARATUS USED:

Diagram



METHOD:

- 1) *Use the pipette to measure 20 ml of Sodium Hydroxide solution into the conical flask.*
- 2) *Add 4 drops of Litmus indicator solution to the conical flask.*
- 3) *Fill the burette with hydrochloric acid and slowly add the acid to the alkali.*
- 4) *Stop adding acid when the solution in the flask just begins to turn red. Record the volume of acid.*
- 5) *Repeat the experiment, and then average the two volumes of acid needed for neutralisation.*

RESULTS:			
Titration number	1	2	Average value
Volume of acid added (ml)			
CONCLUSIONS:			

Modified from Students Laboratory Notebook.